

Clinical Competency Form for Mask Transport



This document is a Clinical Competency Statement: The participant must perform setup and operating proficiency without assistance and/or direct supervision regarding the ANAPOD™ Humi-Therm Heated Humidification System. This document is a self-assessment for registered practitioners.

Applicant Name:			
Hospital/Dept:			
Date:			
	ANAPOD Competency Check Off for Transport	Met	Not Met
1	Establish the need for ANAPOD™		
2	Gather correct equipment: ANAPOD™ Controller ANAPOD™ “Quick Connect” Circuit Cable ANAPOD™ Power Supply Cord ANAPOD™ Humi-Therm Transport Circuit (A9244) Optional: ANAPOD™ Humi-Therm Transport Kit: Adult Oxygen Mask Kit (9439) Pediatric Oxygen Mask Kit (9440) Adult Non-Rebreather Mask Kit (9442) Pediatric Oxygen Mask Kit (9443) Oxygen Tank		
3	Plug power supply cable into back panel of ANAPOD™ and to a power source *Controller DOES NOT have a battery backup power!		
4	Connect Circuit cable to ANAPOD™ (Red and Black Arrows together)		
5	Open transport circuit (A9244) and transport mask kit		
6	Connect 4 pin Circuit Cable to Circuit		
7	Turn on ANAPOD™ and set temperate to desired setting by moving up and down arrows <i>* Set ANAPOD™ 5°C higher than desired patient temp (i.e. To deliver 37°C to patient set ANAPOD™ at 42°C)</i>		
8	Ensure Warming Light comes on <i>*Warming light will remain on until set temperature is reached, then blinks intermittently</i>		
9	Disconnect Circuit Cable – Alarm should sound, AND visual disconnect alarm should be lit		
10	Hit alarm silence – alarm should silence BUT visual alarm should remain lit <i>*Alarm will remain silenced for 2 minutes</i>		
11	Reconnect Circuit Cable – visual alarm should disappear		
12	Inject 40 mL of provided water into circuit using feed port <i>*This ensures equal distribution of water within circuit</i>		
13	Allow circuit to warm up		
14	Confirm oxygen tubing is connected to end of circuit and turn on oxygen flow		
15	Remove Pop-off valve and adapter from circuit		
16	Attach mask, confirm flow through device <i>*Mask transport kit AND transport circuit comes with oxygen tubing, discard circuit oxygen tubing and use transport kit tubing (longer tubing)</i>		

17	Visualize condensation in the circuit, there should be no free-standing water *Condensation must be present, or patient is not receiving humidification		
18	Place device on patient		
19	Continue to monitor patient and condensation within circuit. *Add water to circuit as necessary, do not oversaturate circuit		

Cleaning Instructions			
	Dispose of all single patient items (circuit, adapters, cannula or mask, etc.)		
	Wipe down controller, cable, and bracket (non-disposables) with germicidal wipes *Do not immerse controller in any liquid as it may cause damage.		

I confirm that I have evaluated the above-named person and can state that he/she demonstrates competency in using the named medical device.

Verifier/Manager/Educator Signature:	
Verifier/Manager/Educator Print Name:	

Status: Pass Needs More Training