


# MRI Transport Circuit



**Westmed**  
When Your Care Is Critical

# Features



**Integral PEEP Valve,  
4 – 20 cmH<sub>2</sub>O**



**Also Available with CO<sub>2</sub> Easy™  
Exhaled CO<sub>2</sub> Detector**



**Monitoring Port with  
Tethered Cap**



**MR-Safe**

A complete, lightweight, ready-to-use, non-metallic, non-electrically conductive MRI transport circuit in a single small package. Compatible with Westmed BagEasy™ Manual Resuscitator and with some portable transport ventilators that have a passive exhalation valve. Features a minimal deadspace non-rebreathing valve with a swivel patient connector, adjustable PEEP valve and a proximal monitoring port that can be used for airway pressure or sidestream end-tidal capnography. Now available with Adult CO<sub>2</sub> Easy™ Carbon Dioxide Detector: 2-inch diameter indicator is easy to read and the anti-fog lens reduces condensation from exhaled gas. Purple-to-yellow color change allows exhaled CO<sub>2</sub> to be confirmed from a distance.

**Expandable Tubing.** The packaged length of the 22 mm corrugated tubing is only 28 inches, but the entire circuit expands up to 8 feet to allow convenient ventilator positioning away from the patient or outside of the MRI magnet.

## Specifications.

- PEEP Valve Range (Calibrated Scale) = 4 to 20 cmH<sub>2</sub>O (±4 cmH<sub>2</sub>O); Minimum <3 cmH<sub>2</sub>O; Maximum >20 cmH<sub>2</sub>O.
- Deadspace of non-rebreathing valve = 10 mL.
- Tubing Length = ~28 inches, collapsed; ~8 feet (96 inches), extended.
- Minimum tubing compression/compliance (fully collapsed) = 0.47 mL/cmH<sub>2</sub>O.
- Maximum tubing compression/compliance (fully extended) = 0.95 mL/cmH<sub>2</sub>O.
- Bacterial/viral ventilator outlet filter included.



## Benefits

- Lightweight, decreases pull on ET tube
- MRI-compatible
- Expandable tubing, ~28 – 96 inches
- Minimal Deadspace / Integral PEEP Valve
- Patient swivel connector, easy to attach
- Pre-assembled for simplicity in setup and use

### ORDERING INFORMATION

REF 9317 • MRI Ventilator Circuit, 96" MRI Disposable Circuit with Patient Valve and PEEP, 10/case.

REF 9318 • MRI Ventilator Circuit, 96" MRI Disposable Circuit with Patient Valve, PEEP, CO<sub>2</sub> Easy, 10/case.

Compatible with Westmed BagEasy Manual Resuscitator and with portable transport ventilators that have a passive exhalation valve.

**Westmed**  
When Your Care Is Critical

#### Westmed, Inc.

5580 S. Nogales Highway  
Tucson, AZ 85706

Toll Free: (800) 975-7987  
Phone: (520) 294-7987  
Fax: (520) 294-6061

Distributed by:



Please visit us at  
[www.westmedinc.com](http://www.westmedinc.com)  
for a complete listing of  
Westmed products.